

Anti-Mcu Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Ne86

Item Number: CM0020063



Introduction

Rabbit polyclonal antibody targeting human mitochondrial calcium uniporter MCU. Suitable for WB, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. Supports research on mitochondrial calcium signaling and cellular bioenergetics.

[Learn More](#)

Parameter	Value
Product Name	MCU Rabbit pAb
Remarks/Alias	HsMCU; C10orf42; CCDC109A; MCU
Species	Human
Gene ID (Human)	90550
Gene ID	90550
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 51-170 of human MCU (NP_612366.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NE86
Protein Weight	40kDa
Shipping	Ice bag