

Anti-Mcu Rabbit Polyclonal Antibody For Wb And Elisa - Q8Ne86

Item Number: CM0019927



Introduction

Rabbit polyclonal MCU antibody validated for WB and ELISA. Recognizes human, mouse, rat MCU, a mitochondrial calcium uniporter regulating ATP production, calcium signaling, and cell death pathways. Ideal for mitochondrial research.

[Learn More](#)

Parameter	Value
Product Name	MCU Rabbit pAb
Remarks/Alias	HsMCU; C10orf42; CCDC109A; MCU
Species	Human
Gene ID	90550
Immunogen	Recombinant fusion protein corresponding to amino acids 51-170 of human MCU (NP_612366.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NE86
Protein Weight	40kDa
Shipping	Ice bag