


Anti-Mda5 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Byx4

Item Number: CM0018852



Introduction

Rabbit polyclonal antibody against human MDA5 (IFIH1), a cytoplasmic sensor of viral RNA and activator of antiviral innate immunity. Validated for WB and ELISA applications.

[Learn More](#)

Parameter	Value
Product Name	MDA5 Rabbit pAb
Remarks/Alias	IFIH1; AGS7; Hlcd; IDDM19; MDA-5; MDA5; RLR-2; SGMRT1; interferon induced with helicase C domain 1
Species	Human
GeneID (Human)	64135
GeneID	64135
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9BYX4
Protein Weight	25 kDa/117 kDa
Shipping	Ice bag