


Anti-Mdc1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q14676

Item Number: CM0028405



Introduction

Rabbit polyclonal antibody against human MDC1 (NFB1), a key DNA damage checkpoint protein that binds γ H2AX. Validated in WB, IHC-P, IF/ICC, ELISA. Ideal for DNA repair and cell cycle studies.

[Learn More](#)

Parameter	Value
Product Name	MDC1 Rabbit pAb
Remarks/Alias	NFB1; MDC1
Species	Human
GeneID (Human)	9656
GeneID	9656
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-350 of human MDC1 (NP_055456.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q14676
Protein Weight	227kDa
Shipping	Ice bag