

Anti-Mdh1 Monoclonal Antibody For Wb, Elisa - P40925

Item Number: CM0004680



Introduction

Rabbit monoclonal antibody against human MDH1, validated for Western blot and ELISA. Detects endogenous MDH1 in human, mouse, and rat samples. Suitable for metabolic and cancer research.

[Learn More](#)

Parameter	Value
Product Name	MDH1 Rabbit mAb
Remarks/Alias	KAR; MDHA; MOR2; DEE88; MDH-s; EIEE88; HEL-S-32; MGC:1375; MDH1
Species	Human
GeneID (Human)	4190
GeneID	4190
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 235-334 of human MDH1 (P40925).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40925
Protein Weight	36kDa
Shipping	Ice bag