

Anti-Mdh1 Polyclonal Antibody For Wb And Elisa - P40925

Item Number: CM0028077



Introduction

Recombinant human MDH1 (1-334 aa) polyclonal antibody validated for WB and ELISA in human, mouse, rat. Targets cytoplasmic malate dehydrogenase in malate-aspartate shuttle and TCA cycle.

[Learn More](#)

Parameter	Value
Product Name	MDH1 Rabbit pAb
Remarks/Alias	KAR, MDHA, MOR2, DEE88, MDH-s, EIEE88, HEL-S-32, MGC:1375, MDH1
Species	Human
GeneID (Human)	4190
GeneID	4190
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-334 of human MDH1 (NP_005908.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40925
Protein Weight	36kDa
Shipping	Ice bag