

Anti-Mdh2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P40926

Item Number: CM0003133



Introduction

High-quality rabbit monoclonal antibody targeting human MDH2, validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat MDH2. Ideal for mitochondrial metabolism and neurological disorder research.

[Learn More](#)

Parameter	Value
Product Name	MDH2 Rabbit mAb
Remarks/Alias	MDH; MOR1; DEE51; M-MDH; EIEE51; MGC:3559; MDH2
Species	Human
GeneID (Human)	4191
GeneID	4191
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human MDH2 (P40926)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40926
Protein Weight	36kDa
Shipping	Ice bag