

Anti-Mdm2 Rabbit Monoclonal Antibody For Ihc-P And Elisa - Q00987

Item Number: CM0004384



Introduction

Rabbit monoclonal antibody against human MDM2, validated for IHC-P and ELISA. A synthetic peptide immunogen targets the N-terminal region (aa 101-200). Ideal for studying p53 ubiquitination and oncogenic roles in cancer.

[Learn More](#)

Parameter	Value
Product Name	MDM2 Rabbit mAb
Remarks/Alias	HDMX; LSKB; hdm2; ACTFS; MDM2
Species	Human
Gene ID (Human)	4193
Immunogen	A synthetic peptide corresponding to amino acids 101-200 of human MDM2 (NP_001354919.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q00987
Protein Weight	55kDa
Shipping	Ice bag