

Anti-Mdm2 Rabbit Monoclonal Antibody For Wb, Elisa - Q00987

Item Number: CM0007415



Introduction

Rabbit monoclonal antibody against human MDM2 (HDM2), validated for WB and ELISA. Ideal for studying E3 ubiquitin ligase in p53 pathways, apoptosis, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	MDM2 Rabbit mAb
Remarks/Alias	HDMX; LSKB; hdm2; ACTFS; MDM2
Species	Human
GeneID (Human)	4193
GeneID	4193
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 101-200 of human MDM2 (NP_001354919.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q00987
Protein Weight	11kDa/14kDa/24kDa/30kDa/33kDa/35kDa/48kDa/49kDa/55kDa
Shipping	Ice bag