

Anti-Mdm2 Rabbit Polyclonal Antibody For Wb And Elisa - Q00987

Item Number: CM0016011



Introduction

Rabbit polyclonal antibody against human MDM2 (HDM2), suitable for WB and ELISA. Detects human, mouse, and rat MDM2. Ideal for p53 ubiquitination and cancer research studies. - Q00987

[Learn More](#)

Parameter	Value
Product Name	MDM2 Rabbit pAb
Remarks/Alias	HDMX, LSKB, hdm2, ACTFS, MDM2
Species	Human
GeneID (Human)	4193
GeneID	4193
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-430 of human MDM2 (NP_001354919.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q00987
Protein Weight	55kDa
Shipping	Ice bag