

Anti-Mdm2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q00987

Item Number: CM0017174



Introduction

High-quality rabbit polyclonal antibody against human MDM2, validated for WB, IHC-P, and ELISA. Detects human and mouse MDM2 (55kDa), an E3 ubiquitin ligase and key p53 regulator, ideal for cancer research.

[Learn More](#)

Parameter	Value
Product Name	MDM2 Rabbit pAb
Remarks/Alias	HDMX; LSKB; hdm2; ACTFS; MDM2
Species	Human
Gene ID (Human)	4193
Gene ID	4193
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 231-330 of human MDM2 (NP_001354919.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q00987
Protein Weight	55kDa
Shipping	Ice bag