


Anti-Me2 Polyclonal Antibody For Wb, Elisa - P23368

Item Number: CM0025407



Introduction

Rabbit polyclonal antibody targeting human mitochondrial NAD-dependent malic enzyme (ME2). Validated for Western blot and ELISA, the immunogen corresponds to amino acids 220-479. Ideal for studying mitochondrial metabolism and related diseases.

[Learn More](#)

Parameter	Value
Product Name	ME2 Rabbit pAb
Remarks/Alias	ODS1; ME2
Species	Human
GeneID (Human)	4200
GeneID	4200
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 220-479 of human ME2 (NP_001161807.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P23368
Protein Weight	65kDa
Shipping	Ice bag