


Anti-Me3 (Nadp-Me) Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q16798

Item Number: CM0014388



Introduction

Rabbit polyclonal antibody against human ME3 (NADP-dependent malic enzyme, mitochondrial). Validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Uniprot: Q16798.

[Learn More](#)

Parameter	Value
Product Name	ME3 Rabbit pAb
Remarks/Alias	NADP-ME; ME3
Species	Human
Gene ID (Human)	10873
Gene ID	10873
Immunogen	Recombinant Protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16798
Protein Weight	67kDa
Shipping	Ice bag