

Anti-Meaf6 Polyclonal Antibody For Wb, Elisa - Q9Haf1

Item Number: CM0013106

Introduction

High-quality rabbit polyclonal antibody against MEAF6 (CENP-28, EAF6), suitable for Western blot and ELISA applications. Reacts with human and mouse samples. 22 kDa target.

[Learn More](#)

Parameter	Value
Product Name	MEAF6 Rabbit pAb
Remarks/Alias	EAF6; CENP-28; C1orf149; NY-SAR-91; MEAF6
Species	Human
GeneID (Human)	64769
GeneID	64769
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-191 of human MEAF6 (NP_001257804.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9HAF1
Protein Weight	22kDa
Shipping	Ice bag