


Anti-Meckelin Monoclonal Antibody For Wb, Ihc-P, Elisa - Q5Hya8

Item Number: CM0003857



Introduction

Rabbit monoclonal antibody against human Meckelin (TMEM67), validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for ciliopathy and Wnt signaling research.

[Learn More](#)

Parameter	Value
Product Name	Meckelin Rabbit mAb
Remarks/Alias	MKS3; JBTS6; NPHP11; TNEM67; MECKELIN; Meckelin
Species	Human
GeneID (human)	91147
GeneID	91147
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 800-880 of human Meckelin (NP_714915.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5HYA8
Protein Weight	112kDa
Shipping	Ice bag