


Anti-Mecom Rabbit Polyclonal Antibody For Wb, Elisa - Q03112 / Q13465

Item Number: CM0025582



Introduction

Rabbit polyclonal antibody targeting MECOM (EVI1) for Western blot and ELISA. Validated in human and rat. Suitable for research on transcription regulation, histone methylation, and leukemia.

[Learn More](#)

Parameter	Value
Product Name	MECOM Rabbit pAb
Remarks/Alias	EVI1; MDS1; KMT8E; PRDM3; RUSAT2; MDS1-EVI1; AML1-EVI-1; MECOM
Species	Human
Gene ID (Human)	2122
Gene ID	2122
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human MECOM (NP_004982.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q03112/Q13465
Protein Weight	18kDa/117kDa/118kDa/125kDa/139kDa
Shipping	Ice bag