

Anti-Mecr Polyclonal Antibody For Ihc-P, Elisa - Q9Bv79

Item Number: CM0018946



Introduction

MECR Rabbit pAb (UniProt Q9BV79) validated for IHC-P and ELISA. Reacts with human. Ideal for investigating mitochondrial fatty acid synthesis, iron homeostasis, and MECR-related neurological disorders.

[Learn More](#)

Parameter	Value
Product Name	MECR Rabbit pAb
Remarks/Alias	ETR1; NRB1; CGI-63; FASN2B; DYTOABG; MECR
Species	Human
Gene ID	51102
Immunogen	Recombinant protein corresponding to a sequence within amino acids 54-373 of human MECR (NP_057095.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q9BV79
Protein Weight	40kDa
Shipping	Ice bag