

Anti-Mecr Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Bv79

Item Number: CM0008276



Introduction

Recombinant rabbit monoclonal antibody against human MECR, validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Suitable for mitochondrial fatty acid synthesis and iron homeostasis research.

[Learn More](#)

Parameter	Value
Product Name	MECR Rabbit mAb
Remarks/Alias	ETR1; NRBF1; CGI-63; FASN2B; DYTOABG
Species	Human
GeneID (Human)	51102
GeneID	51102
Immunogen	Recombinant Protein corresponding to a sequence within amino acids 54-373 of human MECR(NP_057095.4).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Bv79
Protein Weight	40kDa
Shipping	Ice bag