

Anti-Med1 Rabbit Polyclonal Antibody For Wb, Elisa - Q15648

Item Number: CM0022941



Introduction

High-quality rabbit polyclonal antibody targeting human MED1 (Mediator complex subunit 1). Validated for WB and ELISA. Suitable for studying RNA polymerase II transcription coactivation and nuclear receptor signaling pathways. Cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	MED1 Rabbit pAb
Remarks/Alias	PBP; CRSP1; RB18A; TRIP2; PPARBP; CRSP200; DRIP205; DRIP230; PPARGBP; TRAP220; MED1
Species	Human
Gene ID (Human)	5469
Gene ID	5469
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human MED1 (NP_004765.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q15648
Protein Weight	168kDa
Shipping	Ice bag