

Anti-Med15 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q96Rn5

Item Number: CM0019672



Introduction

Rabbit polyclonal antibody against human MED15 (ARC105), validated for WB, IF/ICC, ELISA. Detects human, mouse, rat MED15 (~87 kDa). Ideal for transcriptional regulation and Mediator complex research.

[Learn More](#)

Parameter	Value
Product Name	MED15 Rabbit pAb
Remarks/Alias	TIG1; CAG7A; CTG7A; PCQAP; TIG-1; TNRC7; ARC105; MED15
Species	Human
Gene ID (Human)	51586
Gene ID	51586
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MED15 (NP_001003891.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96RN5
Protein Weight	87kDa
Shipping	Ice bag