

Anti-Med17 Polyclonal Antibody For Wb, Elisa - Q9Nvc6

Item Number: CM0029515



Introduction

Rabbit polyclonal antibody against MED17, suitable for WB and ELISA. Validated for human, mouse, and rat samples. Ideal for studying Mediator complex and RNA polymerase II transcription.

[Learn More](#)

Parameter	Value
Product Name	MED17 Rabbit pAb
Remarks/Alias	SRB4; CRSP6; CRSP77; DRIP80; TRAP80; MED17
Species	Human
GeneID (Human)	9440
GeneID	9440
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 250-460 of human MED17 (NP_004259.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NVC6
Protein Weight	73kDa
Shipping	Ice bag