

Anti-Med18 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Bue0

Item Number: CM0014543



Introduction

This rabbit polyclonal antibody against MED18 is validated for WB and ELISA, detecting human and mouse MED18, a Mediator complex subunit essential for RNA polymerase II transcription. Ideal for studying transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	MED18 Rabbit pAb
Remarks/Alias	SRB5; p28b; MED18
Species	Human
GeneID (Human)	54797
GeneID	54797
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BUE0
Protein Weight	24kDa
Shipping	Ice bag