

Anti-Med20 Polyclonal Antibody For Wb, Elisa - Q9H944

Item Number: CM0023537



Introduction

Rabbit polyclonal antibody targeting human MED20, a Mediator complex subunit crucial for RNA polymerase II transcription. Suitable for WB and ELISA; cross-reacts with human, mouse, and rat. Ideal for transcription regulation research.

[Learn More](#)

Parameter	Value
Product Name	MED20 Rabbit pAb
Remarks/Alias	SRB2; TRFP; PRO0213; MED20
Species	Human
GeneID (Human)	9477
GeneID	9477
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-212 of human MED20 (NP_004266.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H944
Protein Weight	23kDa
Shipping	Ice bag