

Anti-Med21 Polyclonal Antibody For Wb, Elisa - Q13503

Item Number: CM0014606



Introduction

Rabbit polyclonal antibody against human MED21 (SRB7), validated for Western blot and ELISA. Mediator complex subunit involved in RNA polymerase II transcription regulation. Suitable for human and mouse samples.

[Learn More](#)

Parameter	Value
Product Name	MED21 Rabbit pAb
Remarks/Alias	SRB7; SURB7; hSrb7; MED21
Species	Human
Gene ID (Human)	9412
Gene ID	9412
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-144 of human MED21 (NP_004255.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13503
Protein Weight	16kDa
Shipping	Ice bag