

Anti-Med24 Polyclonal Antibody For Wb And Elisa - 075448

Item Number: CM0030053



Introduction

MED24 Rabbit Polyclonal Antibody validated for Western Blot (WB) and ELISA. Detects human, mouse, and rat MED24, a 110 kDa Mediator complex subunit. Ideal for transcription regulation research.

[Learn More](#)

Parameter	Value
Product Name	MED24 Rabbit pAb
Remarks/Alias	MED5; CRSP4; ARC100; THRAP4; CRSP100; DRIP100; TRAP100; MED24
Species	Human
GeneID (Human)	9862
GeneID	9862
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human MED24 (NP_055630.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075448
Protein Weight	110kDa
Shipping	Ice bag