

Anti-Med4 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Npj6

Item Number: CM0002266



Introduction

MED4 Rabbit Monoclonal Antibody, validated for WB, IHC-P, ELISA. Recognizes human, mouse, rat MED4, a Mediator complex subunit involved in RNA polymerase II transcription. Ideal for transcription regulation research.

[Learn More](#)

Parameter	Value
Product Name	MED4 Rabbit mAb
Remarks/Alias	ARC36, VDRIP, DRIP36, TRAP36, HSPC126, MED4
Species	Human
GeneID (Human)	29079
GeneID	29079
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 150-250 of human MED4 (Q9NPJ6)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NPJ6
Protein Weight	30kDa
Shipping	Ice bag