

Anti-Med6 Rabbit Polyclonal Antibody For Wb And Elisa - 075586

Item Number: CM0030395



Introduction

MED6 Rabbit Polyclonal Antibody (pAb) validated for WB and ELISA. Detects human Mediator subunit 6 (28kDa). Ideal for transcription research and Mediator complex studies.

[Learn More](#)

Parameter	Value
Product Name	MED6 Rabbit pAb
Remarks/Alias	ARC33; NY-REN-28; MED6
Species	Human
GeneID (Human)	10001
GeneID	10001
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-246 of human MED6 (NP_005457.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	075586
Protein Weight	28kDa
Shipping	Ice bag