

Anti-Med7 Rabbit Polyclonal Antibody For Wb And Elisa - O43513

Item Number: CM0022673



Introduction

MED7 Rabbit Polyclonal Antibody validated for WB and ELISA, recognizing human and mouse MED7. MED7 is a Mediator complex subunit crucial for RNA polymerase II transcription, making it ideal for transcription regulation studies.

[Learn More](#)

Parameter	Value
Product Name	MED7 Rabbit pAb
Remarks/Alias	ARC34; CRSP9; CRSP33; MED7
Species	Human
GeneID (Human)	9443
GeneID	9443
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-233 of human MED7 (NP_004261.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O43513
Protein Weight	27kDa
Shipping	Ice bag