

Anti-Mef2A Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q02078

Item Number: CM0017159



Introduction

Rabbit polyclonal antibody against human MEF2A. Validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat samples. Suitable for muscle development and transcription regulation studies.

[Learn More](#)

Parameter	Value
Product Name	MEF2A Rabbit pAb
Remarks/Alias	mef2; ADCAD1; RSRFC4; RSRFC9; MEF2A
Species	Human
Gene ID (Human)	4205
Gene ID	4205
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MEF2A (NP_001124400.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02078
Protein Weight	55kDa
Shipping	Ice bag