

## Anti-Mef2C Rabbit Polyclonal Antibody For Wb, If, Ihc, Elisa - Q06413

Item Number: CM0016085

### Introduction

MEF2C Rabbit Polyclonal Antibody validated for WB, IF-P, IHC-P, and ELISA. Cross-reacts with Human, Mouse, Rat. Ideal for muscle development, neurobiology, and transcription regulation studies. SWISS: Q06413.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | MEF2C Rabbit pAb                                                                                                   |
| Remarks/Alias    | NEDHSIL; DEL5q14.3; C5DELq14.3; MEF2C                                                                              |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 4208                                                                                                               |
| GeneID           | 4208                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 170-380 of human MEF2C (NP_002388.2) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, IF-P, IHC-P, ELISA                                                                                             |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | Q06413                                                                                                             |
| Protein Weight   | 51kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |