


Anti-Megalin/Lrp2 Monoclonal Antibody For Wb, If, Ihc, Fc And Elisa - P98164

Item Number: CM0006333



Introduction

Recombinant rabbit monoclonal antibody targeting human Megalin/LRP2 C-terminal region. Validated for WB, IF-P, IHC-P, FC (intra), and ELISA. Detects human, mouse, and rat LRP2 (~522 kDa), a multiligand endocytic receptor essential for kidney reabsorption, neural development, and lipoprotein metabolism.

[Learn More](#)

Parameter	Value
Product Name	Megalin/LRP2 Rabbit mAb
Remarks/Alias	DBS; GP330; LRP-2; Megalin/LRP2
Species	Human
GeneID (Human)	4036
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 4446-4655 of human Megalin/LRP2 (NP_004516.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, FC (intra), ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P98164
Protein Weight	522kDa
Shipping	Ice bag