

Anti-Meis1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O00470

Item Number: CM0020600



Introduction

Anti-MEIS1 rabbit polyclonal antibody validated for Western blotting, immunofluorescence/immunocytochemistry, and ELISA. Reacts with human, mouse, and rat MEIS1. Ideal for studying hematopoiesis, transcription regulation, and leukemia.

[Learn More](#)

Parameter	Value
Product Name	MEIS1 Rabbit pAb
Remarks/Alias	MEIS1
Species	Human
Gene ID	4211
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 300-390 of human MEIS1 (NP_002389.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00470
Protein Weight	43kDa
Shipping	Ice bag