

Anti-Meis2 Monoclonal Antibody For Wb, Elisa - O14770

Item Number: CM0008105



Introduction

Rabbit monoclonal antibody targeting human MEIS2, validated for Western blot and ELISA. Detects human and mouse MEIS2. Ideal for transcription regulation and developmental biology studies.

[Learn More](#)

Parameter	Value
Product Name	MEIS2 Rabbit mAb
Remarks/Alias	MRG1; CPCMR; HsT18361; MEIS2
Species	Human
GeneID (human)	4212
GeneID	4212
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MEIS2 (O14770).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O14770
Protein Weight	52kDa
Shipping	Ice bag