

Anti-Meis2 Polyclonal Antibody For Wb, Elisa - O14770

Item Number: CM0026378



Introduction

MEIS2 Rabbit Polyclonal Antibody validated for WB and ELISA with cross-reactivity to mouse and rat. Suitable for investigating transcriptional regulation and developmental biology. Recombinant protein immunogen.

[Learn More](#)

Parameter	Value
Product Name	MEIS2 Rabbit pAb
Remarks/Alias	MRG1; CPCMR; HsT18361; MEIS2
Species	Human
Gene ID	4212
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 351-470 of human MEIS2 (NP_733776.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O14770
Protein Weight	52kDa
Shipping	Ice bag