

Anti-Meis3 Polyclonal Antibody For Western Blot And Elisa - Q99687

Item Number: CM0020173

Introduction

High-quality rabbit polyclonal antibody against human MEIS3, validated for Western blot and ELISA. Suitable for research on transcriptional regulation, pancreatic beta-cell survival, and nuclear localization. Reacts with human, mouse, and rat.

[Learn More](#)

Parameter	Specification
Product Name	MEIS3 Rabbit pAb
Remarks/Alias	MRG2; MEIS3
Species	Human
Gene ID (Human)	56917
Gene ID	56917
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 325-375 of human MEIS3 (NP_001287988.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99687
Protein Weight	41kDa
Shipping	Ice bag