

Anti-Mek1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q02750

Item Number: CM0013680



Introduction

Rabbit polyclonal antibody targeting human MEK1 (MAP2K1) for WB, IHC-P, IF/ICC, and ELISA; cross-reacts with mouse and rat. Designed for MAPK/ERK pathway studies and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	MEK1 Rabbit pAb
Remarks/Alias	MEL; CFC3; MEK1; MKK1; MAPKK1; PRKMK1
Species	Human
Gene ID (Human)	5604
Gene ID	5604
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human MEK1 (NP_002746.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02750
Protein Weight	43kDa
Shipping	Ice bag