

Anti-Mek1 Rabbit Monoclonal Antibody For Wb, If/Icc, Ip, Elisa - Q02750

Item Number: CM0003901



Introduction

Rabbit monoclonal antibody targeting human MEK1 (MAP2K1), a dual-specificity MAP kinase kinase. Suitable for WB, IF/ICC, IP, ELISA; cross-reacts with human, mouse, rat. Ideal for studying the MAPK/ERK signaling pathway (UniProt: Q02750).

[Learn More](#)

Parameter	Value
Product Name	MEK1 Rabbit mAb
Remarks/Alias	MEL; CFC3; MEK1; MKK1; MAPKK1; PRKMK1
Species	Human
GeneID (human)	5604
GeneID	5604
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MEK1 (Q02750).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02750
Protein Weight	43kDa
Shipping	Ice bag