

Anti-Mek1 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - Q02750

Item Number: CM0003100



Introduction

Anti-MEK1 rabbit monoclonal antibody for WB, IP, ELISA in human, mouse, rat. Targets MAP2K1/MEK1 (Q02750), a key MAPK/ERK pathway kinase. Ideal for cancer and signal transduction research.

[Learn More](#)

Parameter	Value
Product Name	MEK1 Rabbit mAb
Remarks/Alias	MEL; CFC3; MEK1; MKK1; MAPKK1; PRKMK1
Species	Human
GeneID (Human)	5604
GeneID	5604
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 150-250 of human MEK1 (Q02750)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02750
Protein Weight	43kDa
Shipping	Ice bag