

Anti-Mek1/Mek2 Monoclonal Antibody For Wb, Elisa - Q02750/P36507

Item Number: CM0009922



Introduction

High-quality rabbit monoclonal antibody targeting human MEK1/MEK2 kinases. Suitable for Western blot and ELISA applications, with cross-reactivity to mouse and rat. Essential tool for MAPK/ERK signaling research.

[Learn More](#)

Parameter	Value
Product Name	MEK1/MEK2 Rabbit mAb
Remarks/Alias	CFC3; MAPKK1; MEK1; MKK1; PRKMK1; MEK1/MEK2
Species	Human
Gene ID (Human)	5604/5605
Gene ID	5604/5605
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-360 of human MEK1/MEK2 (Q02750)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02750/P36507
Protein Weight	44kDa
Shipping	Ice bag