

Anti-Mek1/Mek2 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q02750 / P36507

Item Number: CM0007259



Introduction

This rabbit monoclonal antibody targets MEK1/MEK2, validated for Western blot, IHC-P, and ELISA. It reacts with human, mouse, and rat samples and is ideal for studying the MAPK/ERK signaling pathway.

[Learn More](#)

Parameter	Description
Product Name	MEK1/MEK2 Rabbit mAb
Remarks/Alias	CFC3; MAPKK1; MEK1; MKK1; PRKMK1; MEK1/MEK2
Species	Human
Gene ID (Human)	5604/5605
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 170-270 of human MEK1/MEK2 (NP_002746.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02750 / P36507
Protein Weight	44kDa
Shipping	Ice bag