

Anti-Mek1/Mek2 Rabbit Polyclonal Antibody For Wb, Elisa - Q02750/P36507

Item Number: CM0024325



Introduction

Rabbit polyclonal antibody targeting MEK1/MEK2 (MAP2K1/2), validated for Western blot and ELISA. Cross-reacts with mouse and rat. Ideal for MAPK/ERK signaling pathway studies and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	MEK1/MEK2 Rabbit pAb
Remarks/Alias	MEK1/MEK2
Species	Human
GeneID (human)	5604/5605
GeneID	5604/5605
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 200-300 of human MAP2K1 (NP_002746.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse,Rat
SWISS	Q02750/P36507
Protein Weight	40kDa/43kDa/44kDa
Shipping	Ice bag