

Anti-Mek2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P36507

Item Number: CM0013281



Introduction

Rabbit polyclonal antibody against human MEK2, validated for WB, IHC-P, IF/ICC, and ELISA. Detects MEK2 in human, mouse, and rat. Suitable for MAPK signaling research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] MEK2 Rabbit pAb
Remarks/Alias	CFC4; MEK2; MKK2; MAPKK2; PRKMK2; K2
Species	Human
GeneID (Human)	5605
GeneID	5605
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MEK2 (NP_109587.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P36507
Protein Weight	44kDa
Shipping	Ice bag