

Anti-Mek2 Rabbit Monoclonal Antibody (Ko Validated) For Wb And Elisa - P36507

Item Number: CM0002154



Introduction

Rabbit monoclonal antibody against human MEK2 (MAP2K2), KO validated for Western blotting (WB) and ELISA. Reacts with human and mouse. Predicted molecular weight: 44 kDa. Ideal for studying MAPK signaling.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] MEK2 Rabbit mAb
Remarks/Alias	CFC4; MEK2; MKK2; MAPKK2; PRKMK2; [KO/KD Validated] MEK2
Species	Human
Gene ID (Human)	5605
Gene ID	5605
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MEK2 (P36507)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P36507
Protein Weight	44kDa
Shipping	Ice bag