

Anti-Mek5 Monoclonal Antibody For Wb And Elisa - Q13163

Item Number: CM0002581



Introduction

Anti-MEK5 monoclonal antibody validated for WB and ELISA. Reacts with human, mouse, rat. Targets MAP2K5 kinase involved in MAPK signaling, apoptosis, and cardiac development. Ideal for research use.

[Learn More](#)

Parameter	Value
Product Name	MEK5 Rabbit mAb
Remarks/Alias	MEK5
Species	Human
Gene ID	5607
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 250-350 of human MEK5 (Q13163).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13163
Protein Weight	50kDa
Shipping	Ice bag