

Anti-Mek5 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q13163

Item Number: CM0024843



Introduction

Rabbit polyclonal antibody targeting human MEK5 (MAP2K5), validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for MAPK signaling, apoptosis, and cardiac development research.

[Learn More](#)

Parameter	Value
Product Name	MEK5 Rabbit pAb
Remarks/Alias	MEK5; MAPKK5; PRKMK5; HsT17454
Species	Human
GeneID (Human)	5607
GeneID	5607
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human MEK5 (NP_660143.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13163
Protein Weight	50kDa
Shipping	Ice bag