

Anti-Mek6 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P52564

Item Number: CM0030981



Introduction

MEK6 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA in human, mouse, rat. Recombinant immunogen targeting stress-activated MAPKK6; ideal for cytokine and stress signaling research.

[Learn More](#)

Parameter	Value
Product Name	MEK6 Rabbit pAb
Remarks/Alias	MEK6; MKK6; MAPKK6; PRKMK6; SAPKK3; SAPKK-3
Species	Human
GeneID (Human)	5608
GeneID	5608
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-334 of human MEK6 (NP_002749.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52564
Protein Weight	37kDa
Shipping	Ice bag