


Anti-Mek7 (Map2K7) Monoclonal Antibody For Wb, Elisa - O14733

Item Number: CM0002025



Introduction

Rabbit monoclonal anti-MEK7 antibody validated for WB and ELISA. Reacts with human, mouse, rat. Detects a 47 kDa band. Essential for JNK signaling, stress-induced apoptosis. O14733.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] MEK7 Rabbit mAb
Remarks/Alias	MEK; MKK7; JNKK2; MEK 7; MAPKK7; PRKMK7; SAPKK4; SAPKK-4; K7
Species	Human
GeneID (Human)	5609
GeneID	5609
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 320-419 of human MEK7 (O14733).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14733
Protein Weight	47 kDa
Shipping	Ice bag