


Anti-Mekk2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Y2U5

Item Number: CM0002666



Introduction

High-quality rabbit monoclonal antibody targeting human MEKK2 (MAP3K2), validated for Western blot, IHC-P, and ELISA. Cross-reacts with mouse and rat. Essential tool for MAP kinase research.

[Learn More](#)

Parameter	Value
Product Name	MEKK2 Rabbit mAb
Remarks/Alias	MEKK2; MEKK2B
Species	Human
GenelD (human)	10746
GenelD	10746
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y2U5
Protein Weight	70kDa
Shipping	Ice bag