

Anti-Mekk3 (Map3K3) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q99759

Item Number: CM0019547



Introduction

Rabbit polyclonal antibody against MEKK3 (MAP3K3) for Western blotting, IHC-P, and ELISA. Validated in human, mouse, and rat samples; suitable for studying MAPK signaling and transcription factor activation.

[Learn More](#)

Parameter	Value
Product Name	MEKK3 Rabbit pAb
Remarks/Alias	MEKK3; MAPKKK3
Species	Human
Gene ID	4215
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human MEKK3 (NP_976226.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99759
Protein Weight	71kDa
Shipping	Ice bag