

Anti-Mekk3 Monoclonal Antibody For Wb - Q99759

Item Number: CM0012719



Introduction

Rabbit monoclonal antibody against human MEKK3 (MAP3K3), validated for WB and ELISA, reactive with human and mouse, protein weight 71kDa. Ideal for signal transduction research.

[Learn More](#)

Parameter	Value
Product Name	MEKK3 Rabbit mAb
Remarks/Alias	MEKK3; MAPKKK3
Species	Human
Gene ID (Human)	4215
Gene ID	4215
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MEKK3 (Q99759).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q99759
Protein Weight	71kDa
Shipping	Ice bag